

Calculation:

The experiment is divided into two parts.

1. Determination of end point or equivalent point
2. Determination of pKa. (Using table-1 corresponding to pKa find out which weak acid is present)

Table-1

Acid	pKa	Gram equivalent weight
CH ₃ COOH (acetic)	4.7	60/1 = 60
HCOOH (formic)	3.7	46/1 = 46
H ₂ CO ₃ (carbonic)	6.3	62/2 = 31
H ₂ C ₂ O ₄ (oxalic)	1.2	126/2 = 63

From the equivalence point find out the normality of weak acid using the formula

$$(NV)_{\text{acid}} = (NV)_{\text{NaOH}}$$


$$N_{\text{acid}} = \frac{(N'V')_{\text{NaOH}}}{V_{\text{acid}}}$$

Amount of acid present in 1000 cm³ of its solution = N_{acid} × gram equivalent weight of acid


Model Procedure /Flow Chart:

- 1) Take 50 cm³ of soft drink which in this case is CH₃COOH in a clean 250 cm³ beaker
- 2) Insert the glass calomel electrode such that it touches the bottom and connect it to pH meter
- 3) Read the initial pH and note down
- 4) Fill the burette with sodium hydroxide
- 5) Note the initial reading and keep adding 0.2 cm³ of NaOH at a time and take reading from pH meter
- 6) Continue the process till there is sudden jump in pH
- 7) Take five more readings and plot the graph.

Model graph-1:



Model graph-2:



Model Calculation:

For equivalence point,

Normality of a weak acid is

$$(NV)_{\text{acid}} = (NV)_{\text{NaOH}}$$

$$N_{\text{acid}} = \frac{(NV)_{\text{NaOH}}}{V_{\text{acid}}}$$

Amount of acid present in 1000cm^3 of its solution

$$= N_{\text{acid}} \times \text{Gram molecular weight of acid}$$


Tabulation:


Volume of NaOH added in cm ³	pH	ΔpH	ΔV	$\frac{\Delta pH}{\Delta V}$
0.0	3.37	-	-	-
0.2	3.99	0.62	0.2	3.1
0.4	4.36	0.37	0.2	1.85
0.6	4.59	0.23	0.2	1.15
0.8	4.86	0.27	0.2	1.35
1.0	5.12	0.26	0.2	1.30
1.2	5.48	0.36	0.2	1.80
1.4	6.28 6.18	0.80	0.2	4.4
1.6	10.35 10.35	4.07	0.2	20.35
1.8	11.35	1.00	0.2	5
2.0	11.61	0.26	0.2	1.3
2.2	11.73	0.12	0.2	0.6
2.4	11.83	0.10	0.2	0.5
2.6	11.92	0.09	0.2	0.45
2.8	11.99	0.07	0.2	0.35
3.0	12.04	0.05	0.2	0.25
3.2	12.09	0.05	0.2	0.25
3.4	12.14 12.13	0.05	0.2	0.25
3.6	12.20	0.06	0.2	0.30

Scale

$$x\text{-axis} \rightarrow 1 \text{ unit} = 0.2 \text{ cm}^3$$

$$y\text{-axis} \rightarrow 1 \text{ unit} = 1 \text{ M}$$





Calculation:

$$(NV)_{\text{acid}} = (NV)_{\text{NaOH}}$$

$$N_{\text{acid}} = \frac{(NV)_{\text{NaOH}}}{V_{\text{acid}}} = \frac{0.05 \times 1.8}{50} = 0.0018 \text{ Molar}$$

Amount of acid present in 100 cm^3 of solution

$$\begin{aligned} &= N_{\text{acid}} \times \text{Gram equivalent weight of acid} \\ &= (0.0018 \times 60) \text{ gm} = 0.108 \text{ gm} \approx \boxed{0.11 \text{ gm}} \end{aligned}$$

Inference:

In this experiment we can measure the pH of solution which lies between 2.5 to 3.5 which is very good. By doing this experiment [using a combined glass calomel electrode formed due to addition of a strong base & weak acid] In endpoint, acid is neutralized.

Relevance to Society & Environment:

We should be aware of pH of water. ~~of~~ Very low pH can cause damage to our teeth and can also have various side effects. High content of acid in softdrinks and water is harmful.

Report: Total acidity of the given soft drink = 0.11 gm
Acid content

Evaluation of Experiment - 3		
Components	Marks	
	Max	Obtained
Model Procedure, Model Graph & Calculation	16	16
Equivalence Point & Execution	20	20
Inference & Societal Relevance	04	03
Total	40	39

Signature of Teacher

[Signature]